



United 667

MIGHTY DEFENSE Disinfecting Wipes

BACTERICIDAL-TUBERCULOCIDAL-
FUNGICIDAL-VIRUCIDAL*

Kills Pandemic 2009 H1N1 Influenza A virus

Disinfects in 1 minute.

ACTIVE INGREDIENTS:

Ethyl Alcohol 72.50%

Didecyl Dimethyl Ammonium Chloride..... 0.33%

OTHER INGREDIENTS*: 27.17%

TOTAL: 100.00%

*Does not include the weight of the dry wipes.

KEEP OUT OF REACH OF CHILDREN CAUTION

See side panels for additional
precautionary statements.

EPA Reg. No. 88494-4-9250 EPA Est. No. 9250-IL-01

CASE NET CONTENTS: 6 x 1 Canisters

NET CONTENTS PER CANNISTER: 1 lb. 3.5 ozs.

150 Pre-Saturated Wipes

6"x 6.75" (15.24cm x 17.15cm)

Manufactured and Sold By:
UNITED LABORATORIES, INC.
320 37th Avenue • St. Charles, IL 60174
1-800-323-2594 • www.unitedlabsinc.com

DO NOT FLUSH



EFFICACY TESTS HAVE DEMONSTRATED THAT THIS PRODUCT IS AN EFFECTIVE, VIRUCIDE, BACTERICIDE, TUBERCULOCIDE AND FUNGICIDE IN THE PRESENCE OF [5%] ORGANIC SOIL.

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

This product is not for use on critical and semi-critical medical device surfaces.

CLEANING: Use one or more wipes, as necessary, to wet surface sufficiently and to thoroughly clean the surface.

DISINFECTING: To disinfect hard, non-porous surfaces, use one or more wipes, as necessary, to thoroughly wet the surface to be treated. Treated surface must remain visibly wet for one minute to achieve complete disinfection [or kill] of all pathogens listed on this label. For visibly soiled surfaces, clean before following disinfecting instructions. Allow surfaces to air dry.

United 667 MIGHTY DEFENSE Disinfecting Wipes are intended for use on hard non-porous surfaces in: Hospitals and other institutional acute and long term Healthcare facilities, Ambulatory Care Centers, Medical/Dental/Veterinary Offices, Airports, Athletic Facilities, Gyms, Health Clubs, Schools, Transportation Terminals, Vehicles, Aircraft, Schools, Childcare Facilities, Daycare Centers, Nurseries, Libraries, Business offices and manufacturing facilities, Police/Fire/EMS facilities and vehicles, Government offices, Military bases, facilities, vehicles, and aircraft, Detention and Correctional Facilities and vehicles, Animal Housing/Holding Facilities, Animal Hospitals, Animal Life Science Laboratories, Animal Transportation Vehicles, Breeding Establishments, Grooming Establishments, Kennels, Pet Shops, Tack Shops, Veterinary Clinics, Veterinary Examination Rooms, Veterinary Operating Rooms, Veterinary Waiting Rooms, Veterinary X-Ray Rooms, Zoos.

***VIRUCIDAL EFFICACY AGAINST THE FOLLOWING ORGANISMS IN ONE MINUTE:** This product is an effective virucide on hard, non-porous surfaces, in the presence of an organic load, against Adenovirus type 5, Avian Influenza A Virus, Cytomegalovirus, Enterovirus EV-D68, Hepatitis B Virus, Hepatitis C Virus, Herpes Simplex Virus Type 1, Human Coronavirus, Human Immunodeficiency Virus (HIV-1), Influenza Virus type A₂, Influenza B, Norovirus, Poliovirus serotype 1, Respiratory Syncytial Virus, Rhinovirus, Rotavirus Strain WA, and Vaccinia virus, when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry.

KILLS HIV-1, HBV and HCV ON PRE-CLEANED ENVIRONMENTAL SURFACES/OBJECTS PREVIOUSLY SOILED WITH BLOOD/BODY FLUIDS in health care settings Hospitals, Nursing Homes or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1, Hepatitis B Virus, and Hepatitis C Virus.

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV, and HCV FOR SURFACES OR OBJECTS SOILED WITH BLOOD OR BODY FLUIDS.

PERSONAL PROTECTION: When handling items soiled with blood or body fluids, use disposable latex gloves, gowns, masks and eye coverings.

CLEANING PROCEDURES: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of **United 667 MIGHTY DEFENSE Disinfecting Wipes**.

CONTACT TIME: Allow surface to remain wet for 1 full minute.

DISPOSAL OF INFECTIOUS MATERIALS: Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious waste disposal.

BACTERICIDAL AGAINST THE FOLLOWING ORGANISMS IN ONE MINUTE: DRUG RESISTANT AND MULTI-DRUG RESISTANT BACTERIA: **United 667 MIGHTY DEFENSE Disinfecting Wipes** are an effective disinfectant on hard non-porous surfaces, in the presence of an organic load, against *Acinetobacter baumannii* – Multi-Drug resistant, *Carbapenem resistant Enterobacter cloacae* (CRE) (KPC), *Carbapenem resistant Enterobacter cloacae* (CRE) (New Delhi), *Enterococcus faecium* – Vancomycin resistant (VRE), *Escherichia coli* – ESBL Producing, *Klebsiella pneumoniae* – ESBL Producing, *Staphylococcus aureus* – Methicillin resistant (MRSA), *Staphylococcus aureus* – Vancomycin resistant (VRSA), *Staphylococcus aureus* Panel – Vancomycin intermediate (VISA), when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry.

BACTERIA: **United 667 MIGHTY DEFENSE Disinfecting Wipes** are an effective disinfectant on hard non-porous surfaces, in the presence of an organic load, against *Bordetella pertussis*, *Burkholderia cepacia*, *Campylobacter jejuni*, *Enterobacter aerogenes*, *Escherichia coli* (E. coli) O157:H7, *Klebsiella pneumoniae*, *Legionella pneumophila*, *Listeria monocytogenes*, *Proteus vulgaris*, *Pseudomonas aeruginosa*, *Salmonella enterica* (Salmonella), *Shigella dysenteriae*, *Staphylococcus aureus*, *Streptococcus pyogenes* (Strep), and *Streptococcus pneumoniae*, when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry.

TUBERCULOCIDAL: **United 667 MIGHTY DEFENSE Disinfecting Wipes** are an effective disinfectant on hard non-porous surfaces, in the presence of an organic load, against *Mycobacterium bovis*, when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry.

FUNGICIDAL: **United 667 MIGHTY DEFENSE Disinfecting Wipes** are an effective fungicide on hard non-porous surfaces against *Trichophyton interdigitale* (causative agent of Athlete's Foot Fungus), when used as directed, and *Candida albicans*, when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry.

ANIMAL PATHOGENS: **United 667 MIGHTY DEFENSE Disinfecting Wipes** are an effective disinfectant on hard non-porous surfaces against Avian Influenza A Virus, *Bordetella bronchiseptica*, Bovine Diarrhea Virus, Canine distemper Virus, Canine parvovirus, Duck Hepatitis B Virus, Feline calicivirus, and Feline leukemia Virus.

Veterinary Clinics, Animal Life Science Laboratories, Pet Shops: Remove all animals and feed from premises, vehicles, and enclosures. Remove all litter and manure from floors, walls, and surfaces of facilities and fixtures occupied or traversed by animals. Empty all troughs, racks and other feeding and watering appliances. Thoroughly clean all surfaces with soap or detergent and rinse with water. Use one or more wipes to wet surfaces for 1 minute. Immerse all handling and restraining equipment as well as forks, shovels, and scrapers used for removing litter and manure. Ventilate buildings, vehicles, and other closed spaces. Do not house animals or employ equipment until treatment has been absorbed, set, or dried. Thoroughly scrub treated feed troughs, racks and other feeding and watering appliances with soap or detergent, then rinse with potable water before reuse.

DISPENSING INSTRUCTIONS

To Dispense Wipes: Open lid, remove entire cover, remove and discard seal. From center of wipes roll pull up a corner of the wipe, twist it and thread it through the hole in the container cover. Pull through about one inch. Replace cover tightly. Pull out first wipe and tear off at a 45° angle. Remaining wipes feed automatically, ready for next use. When not in use, keep the tip of the exposed wipe inside the well in the cover, and the lid snapped closed.

STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage or disposal.

PESTICIDE STORAGE: Always store this product and all pesticides in the original container. If a leaky container must be contained within another, mark the outer container to identify the contents. Protect from extreme heat or cold.

PESTICIDE DISPOSAL: Waste resulting from the use of this product must be disposed of on site or at an approved waste disposal facility.

CONTAINER DISPOSAL: Nonrefillable container. Do not reuse or refill this container. Offer for recycling if available, or discard empty container in trash.

PRECAUTIONARY STATEMENTS

Hazards to Humans and Domestic Animals:

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using the toilet.

FIRST AID

IF IN EYES: Hold eyes open and rinse slowly and gently with water for 15–20 minutes. Remove contact lenses if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice.

Have the product container or label with you when calling a poison control center, doctor, or going for treatment. For emergency information concerning this product, call the National Pesticides Information Center [NPIC] at 1-800-858-7378 seven days a week, 6:30 am to 4:30 pm Pacific Time [NPIC Web site: www.npic.orst.edu].

PHYSICAL OR CHEMICAL HAZARDS:

Flammable. Keep away from heat or open flame.

United 667 MIGHTY DEFENSE Disinfecting Wipes will not harm most surfaces: including stainless steel, aluminum, other hard and soft metals, sealed quartz or other solid man-made counter top materials, most plastics, fiberglass, vinyl, nylon, and other non-porous upholstery materials used in a medical or veterinary or dental environment, and most painted surfaces. For plastic and painted surfaces, spot test on an inconspicuous area before use. Not recommended for use on acrylic plastics or coated optical surfaces. Do not use on unpainted wood, or leather. Not for use on people or pets.

No warranty, either expressed or implied, is given to the usage contrary to the instructions given herein.

Made in USA

625



DO NOT FLUSH



United 667

MIGHTY DEFENSE

Disinfecting Wipes

BACTERICIDAL-TUBERCULOCIDAL-FUNGICIDAL-VIRUCIDAL*

Kills Pandemic 2009 H1N1 Influenza A virus
Disinfects in 1 minute.

ACTIVE INGREDIENTS:

| | |
|---|---------|
| Ethyl Alcohol..... | 72.50% |
| Didecyl Dimethyl Ammonium Chloride..... | 0.33% |
| OTHER INGREDIENTS**: | 27.17% |
| TOTAL: | 100.00% |

**Does not include the weight of the dry wipes.

KEEP OUT OF REACH OF CHILDREN CAUTION

See side panels for additional precautionary statements.

EPA Reg. No. 88494-4-9250 EPA Est. No. 9250-IL-01

NET CONTENTS: 1 lb. 3.5 ozs.

150 Pre-Saturated Wipes 6"x 6.75" (15.24cm x 17.15cm)

Manufactured and Sold By:

UNITED LABORATORIES, INC.

320 37th Avenue • St. Charles, IL 60174

1-800-323-2594 • www.unitedlabsinc.com

EFFICACY TESTS HAVE DEMONSTRATED THAT THIS PRODUCT IS AN EFFECTIVE, VIRUCIDE, BACTERICIDE, TUBERCULOCIDE AND FUNGICIDE IN THE PRESENCE OF [5%] ORGANIC SOIL.

DIRECTIONS FOR USE

It is a violation of Federal Law to use this product in a manner inconsistent with its labeling.

This product is not for use on critical and semi-critical medical device surfaces.

CLEANING: Use one or more wipes, as necessary, to wet surface sufficiently and to thoroughly clean the surface.

DISINFECTING: To disinfect hard, non-porous surfaces, use one or more wipes, as necessary, to thoroughly wet the surface to be treated. Treated surface must remain visibly wet for one minute to achieve complete disinfection [or kill] of all pathogens listed on this label. For visibly soiled surfaces, clean before following disinfecting instructions. Allow surfaces to air dry.

United 667 MIGHTY DEFENSE Disinfecting Wipes are intended for use on hard non-porous surfaces in: Hospitals and other institutional acute and long term Healthcare facilities, Ambulatory Care Centers, Medical/Dental/Veterinary Offices. Airports, Athletic Facilities, Gyms, Health Clubs, Schools, Transportation Terminals, Vehicles, Aircraft. Schools, Childcare Facilities, Daycare Centers, Nurseries. Libraries. Business offices and manufacturing facilities. Police/Fire/EMS facilities and vehicles. Government offices. Military bases, facilities, vehicles, and aircraft. Detention and Correctional Facilities and vehicles. Animal Housing/Holding Facilities. Animal Hospitals, Animal Life Science Laboratories, Animal Transportation Vehicles, Breeding Establishments, Grooming Establishments, Kennels, Pet Shops, Tack Shops, Veterinary Clinics, Veterinary Examination Rooms, Veterinary Operating Rooms, Veterinary Waiting Rooms, Veterinary X-Ray Rooms, Zoos.

***VIRUCIDAL EFFICACY AGAINST THE FOLLOWING ORGANISMS IN ONE MINUTE:** This product is an effective virucide on hard, non-porous surfaces, in the presence of an organic load, against *Adenovirus type 5*, *Avian Influenza A Virus*, *Cytomegalovirus*, *Enterovirus EV-D68*, *Hepatitis B Virus*, *Hepatitis C Virus*, *Herpes Simplex Virus Type 1*, *Human Coronavirus*, *Human Immunodeficiency Virus (HIV-1)*, *Influenza Virus type A₂*, *Influenza B*, *Norovirus*, *Poliovirus serotype 1*, *Respiratory Syncytial Virus*, *Rhinovirus*, *Rotavirus Strain WA*, and *Vaccinia virus*, when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry.

KILLS HIV-1, HBV and HCV ON PRE-CLEANED ENVIRONMENTAL SURFACES/ OBJECTS PREVIOUSLY SOILED WITH BLOOD/ BODY FLUIDS in health care settings Hospitals, Nursing Homes or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood or body fluids, and in which the surfaces/objects likely to be soiled with blood or body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1, Hepatitis B Virus, and Hepatitis C Virus.

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV, and HCV FOR SURFACES OR OBJECTS SOILED WITH BLOOD OR BODY FLUIDS.

PERSONAL PROTECTION: When handling items soiled with blood or body fluids, use disposable latex gloves, gowns, masks and eye coverings.

CLEANING PROCEDURES: Blood and other body fluids must be thoroughly cleaned from surfaces and objects before application of **United 667 MIGHTY DEFENSE Disinfecting Wipes**.

CONTACT TIME: Allow surface to remain wet for 1 full minute.

DISPOSAL OF INFECTIOUS MATERIALS: Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious waste disposal.

BACTERICIDAL AGAINST THE FOLLOWING ORGANISMS IN ONE MINUTE: DRUG RESISTANT AND MULTI-DRUG RESISTANT BACTERIA: **United 667 MIGHTY DEFENSE Disinfecting Wipes** are an effective disinfectant on hard non-porous surfaces, in the presence of an organic load, against *Acinetobacter baumannii* – Multi-Drug resistant, *Carbapenem resistant Enterobacter cloacae* (CRE) (KPC), *Carbapenem resistant Enterobacter cloacae* (CRE) (New Delhi), *Enterococcus faecium* – Vancomycin resistant (VRE), *Escherichia coli* – ESBL Producing, *Klebsiella pneumoniae* – ESBL Producing, *Staphylococcus aureus* – Methicillin resistant (MRSA), *Staphylococcus aureus* – Vancomycin resistant (VRSA), *Staphylococcus aureus* Panel – Vancomycin intermediate (VISA), when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry.

BACTERIA: **United 667 MIGHTY DEFENSE Disinfecting Wipes** are an effective disinfectant on hard non-porous surfaces, in the presence of an organic load, against *Bordetella pertussis*, *Burkholderia cepacia*, *Campylobacter jejuni*, *Enterobacter aerogenes*, *Escherichia coli* (E. coli) O157:H7, *Klebsiella pneumoniae*, *Legionella pneumophila*, *Listeria*

monocytogenes, *Proteus vulgaris*, *Pseudomonas aeruginosa*, *Salmonella enterica* (Salmonella), *Shigella dysenteriae*, *Staphylococcus aureus*, *Streptococcus pyogenes* (Strep), and *Streptococcus pneumoniae*, when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry.

TUBERCULOCIDAL: **United 667 MIGHTY DEFENSE Disinfecting Wipes** are an effective disinfectant on hard non-porous surfaces, in the presence of an organic load, against *Mycobacterium bovis*, when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry.

FUNGICIDAL: **United 667 MIGHTY DEFENSE Disinfecting Wipes** are an effective fungicide on hard non-porous surfaces against *Trichophyton interdigitale* (causative agent of Athlete's Foot Fungus), when used as directed, and *Candida albicans*, when the treated surface is allowed to remain wet for 1 minute. Allow surfaces to air dry.

ANIMAL PATHOGENS: **United 667 MIGHTY DEFENSE Disinfecting Wipes** are an effective disinfectant on hard non-porous surfaces against *Avian Influenza A Virus*, *Bordetella bronchiseptica*, *Bovine Diarrhea Virus*, *Canine distemper Virus*, *Canine parvovirus*, *Duck Hepatitis B Virus*, *Feline calicivirus*, and *Feline leukemia Virus*.

Veterinary Clinics, Animal Life Science Laboratories, Pet Shops: Remove all animals and feed from premises, vehicles, and enclosures. Remove all litter and manure from floors, walls, and surfaces of facilities and fixtures occupied or traversed by animals. Empty all troughs, racks and other feeding and watering appliances. Thoroughly clean all surfaces with soap or detergent and rinse with water. Use one or more wipes to wet surfaces for 1 minute. Immerse all handling and restraining equipment as well as forks, shovels, and scrapers used for removing litter and manure. Ventilate buildings, vehicles, and other closed spaces. Do not house animals or employ equipment until treatment has been absorbed, set, or dried. Thoroughly scrub treated feed troughs, racks and other feeding and watering appliances with soap or detergent, then rinse with potable water before reuse.

DISPENSING INSTRUCTIONS

To Dispense Wipes: Open lid, remove entire cover, remove and discard seal. From center of wipes roll pull up a corner of the wipe, twist it and thread it through the hole in the container cover. Pull through about one inch. Replace cover tightly. Pull out first wipe and tear off at a 45° angle. Remaining wipes feed automatically, ready for next use. When not in use, keep the tip of the exposed wipe inside the well in the cover, and the lid snapped closed.

STORAGE AND DISPOSAL

Do not contaminate water, food or feed by storage or disposal.

PESTICIDE STORAGE: Always store this product and all pesticides in the original container. If a leaky container must be contained within another, mark the outer container to identify the contents. Protect from extreme heat or cold.

PESTICIDE DISPOSAL: Waste resulting from the use of this product must be disposed of on site or at an approved waste disposal facility.

CONTAINER DISPOSAL: Nonrefillable container. Do not reuse or refill this container. Offer for recycling if available, or discard empty container in trash.

PRECAUTIONARY STATEMENTS

Hazards to Humans and Domestic Animals:

CAUTION: Causes moderate eye irritation. Avoid contact with eyes or clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco, or using the toilet.

FIRST AID

IF IN EYES: Hold eyes open and rinse slowly and gently with water for 15–20 minutes. Remove contact lenses if present, after the first 5 minutes, then continue rinsing. Call a poison control center or doctor for treatment advice.

Have the product container or label with you when calling a poison control center, doctor, or going for treatment. For emergency information concerning this product, call the National Pesticides Information Center [NPIC] at 1-800-858-7378 seven days a week, 6:30 am to 4:30 pm Pacific Time [NPIC Web site: www.npic.orst.edu].

PHYSICAL OR CHEMICAL HAZARDS:

Flammable. Keep away from heat or open flame.

United 667 MIGHTY DEFENSE Disinfecting Wipes will not harm most surfaces: including stainless steel, aluminum, other hard and soft metals, sealed quartz or other solid man-made counter top materials, most plastics, fiberglass, vinyl, nylon, and other non-porous upholstery materials used in a medical or veterinary or dental environment, and most painted surfaces. For plastic and painted surfaces, spot test on an inconspicuous area before use. Not recommended for use on acrylic plastics or coated optical surfaces. Do not use on unpainted wood, or leather. Not for use on people or pets.

No warranty, either expressed or implied, is given to the usage contrary to the instructions given herein.

Made in USA

625